

## Hypersensitive Reverse Transcriptase Instructions

### Composition

Cat. No.	8106050	8106200
Hypersensitive Reverse Transcriptase	50 $\mu$ l	200 $\mu$ l
5 $\times$ RT Buffer	0.5 ml	1 ml
RNase-free Water	1.5 ml	1.5 ml
Instructions	1	1

### Storage

It can be stored for up to two years at -20°C.

### Technical Support

R&D Department, Hangzhou Simgen Biotechnology Co., Ltd. E-mail: technical@simgen.cn, Tel: 400-0099-857.

### Introduction

Hypersensitive Reverse Transcriptase is the third generation M-MLV reverse transcriptase obtained by gene modification and recombination technology. Compared with wild-type M-MLV reverse enzyme, this enzyme removes RNase H activity and significantly improves reverse transcription speed and thermal stability (maximum tolerance temperature is 60°C), and enhances tolerance to the complex secondary structure of RNA. Highly sensitive reverse transcriptase has strong anti-interference and very high reverse transcriptase efficiency for low-copy RNA template, and is the preferred reverse transcriptase for RT-PCR nucleic acid detection kit.

### Applications

1. cDNA library construction.
2. RT-PCR and Real Time RT-PCR.
3. Primer extension.
4. RNA sequencing.
5. One-step RT-PCR

### Unit Definition

Product concentration of 200 U/ $\mu$ l. Definition of activity unit: Using Poly (A) as template and Oligo (dT) as primer, the amount of enzyme required to catalyze the incorporation of 1 nmol dTTP within 10 min at 37°C was defined as 1 activity unit (U).

### Purity

The purity of the product is more than 90% by Coomassie brilliant blue staining SDS-PAGE. The product is free of endonuclease, exonuclease and RNase contamination.

### Equipment and Reagents to Be Supplied by User

1. oligo(dT)<sub>12-18</sub> (25  $\mu$ M) or random primers (25  $\mu$ M) or gene specific primers (1  $\mu$ M).
2. dNTPs (10 mM each, Simgen Cat.No.7701100).
3. RNase Inhibitor (Simgen Cat.No.8008125) may be required.
4. RNase-free microcentrifuge tubes .
5. Pipettes and tips (To avoid RNase contamination, RNase-free tips with filter must be selected).
6. Protective equipment such as disposable latex gloves and paper towels.
7. Water bath.
8. A laboratory without RNase: Because RNase is present in saliva and skin, please wear latex gloves and a mask during the whole experiments.

## Protocol

### 1. Add the following reagents to a RNase-free sterilized microcentrifuge tube.

- 1) 2  $\mu$ l oligo(dT)<sub>12-18</sub>(25  $\mu$ M) or 2  $\mu$ l random primers (25  $\mu$ M) or 2  $\mu$ l gene specific primers (1  $\mu$ M).
- 2) 0.5-5  $\mu$ g total RNA or 50-500 ng mRNA;

\* When total RNA is less than 0.5  $\mu$ g, it was recommended to add 1  $\mu$ l RNase Inhibitor (Cat. No. 8008125).

\* If the RNA template needs to be incubated at 70°C for 5 min to destroy the secondary structure, the addition of RNase Inhibitor should not be omitted.

- 3) 1  $\mu$ l dNTPs (10 mM each).
- 4) Add RNase-free Water to 15  $\mu$ l.

\* If the RNA template is GC-rich or has a complex secondary structure, add the following steps: Incubated at 70°C for 5 min to destroy the RNA secondary structure, then quickly place on ice to prevent the secondary structure from re-forming, then briefly centrifuge.

### 2. Add the reagents according to the table below:

The liquid mixture in step 1	15 $\mu$ l
5 $\times$ RT Buffer	4 $\mu$ l
Hypersensitive M-MLV Reverse Transcriptase	1 $\mu$ l *
Total	20 $\mu$ l

\* When total RNA is less than 0.5  $\mu$ g (such as viral RNA reverse transcription), the amount of Hypersensitive M-MLV Reverse Transcriptase should be reduced to 0.05~0.5  $\mu$ l, otherwise it may lead to subsequent PCR amplification to produce non-specific amplification products.

3. Mix gently, if using random primers as primers, hold at 25°C for 10 min.
4. Incubate at 50°C for 30 min.
5. Incubate at 95°C for 5 min, then cool on ice or store below -20°C for later use.
6. Dilute to 50  $\mu$ l with RNase-free Water and take 2~5  $\mu$ l for PCR amplification.